Phone: 402-717-5227 Fax: 402-717-5252



Anatomic Pathology Specimen Identification and Labeling

Cytology (Non-Gynecological) Specimens

- 1. All **specimen containers** for fluids, urine, brushings, etc. <u>MUST</u> be labeled with the following:
 - a. Patient's full legal first and last names.
 - b. Date of birth OR
 - c. Clinic medical record number.
 - d. Site/source of fluid, including laterality (right or left) as applicable.
 - 1) Note: Do <u>NOT</u> use name of procedure as site/source (e.g. use "abdominal" <u>NOT</u> "paracentesis" fluid; use "pleural" <u>NOT</u> "thoracentesis" fluid, etc.).
 - e. Date and time of collection.
- 2. Additionally, any container <u>MUST</u> be initialed by collecting personnel.
 - a. Specimens with removable lids need the primary container labeled, not the lid.
- 3. Any **glass slide** is considered to be a "primary container", and <u>requires</u> the following:
 - a. Patient's full legal first and last names written in **pencil** on the frosted end of slide.
 - 1) Note: Do <u>NOT</u> use ink or marker (washes off during processing).
 - b. Date of birth OR
 - c. Clinic medical record number.
 - d. Site/source, including laterality (right or left) and laterality (upper, lower, inferior, superior, etc.) as applicable.
- 4. Slides must also indicate whether they were air dried or spray fixed.
 - a. Write "AD" for air dried and "Fix" for fixed as applicable on slides under name, date of birth/MR#, and site/source.
- 5. A patient's registration label may be placed on the outside of a slide holder in addition to labeling individual slides.
 - a. It is <u>NOT</u> acceptable to label the slide holder in lieu of labeling individual slides.
 - b. Unlabeled slides will **NOT** be accepted.
- 6. All **test requisitions** MUST contain the following:
 - a. Patient's full legal first and last names.
 - b. Date of birth.
 - c. Clinic medical record number.
 - d. Site/source of specimen.
 - e. Date and time of collection.
 - f. First initial and full last name of collecting personnel.
- 7. Information on specimen containers/slides must exactly match information on test requisitions.
- 8. Identification and/or labeling issues may result in delay of testing and results.
- 9. Refer to the Cytology Collection Guide on CHI Health Laboratory website for specific information on specimen requirements and collection techniques.